

# FELINE LEUKEMIA VIRUS ANTIGEN & FELINE IMMUNODEFICIENCY ANTIBODY TEST KIT

## RapidSTATUS™ FIV/FeLV

FOR EXPORT ONLY

### TEST PRINCIPLE

The RapidSTATUS™ FIV/FeLV test is a rapid status test for detection of feline immunodeficiency virus antibody and leukemia virus antigen in serum, plasma or anticoagulated whole blood.

### KIT CONTENTS

- Test Devices
- Pipettes
- Chase Buffer
- Instructions for use

### STORAGE

- The kit should be stored at 2-30°C.
- Avoid freezing.
- Avoid storage in direct sunlight.

### PRECAUTIONS

- For veterinary use only.
- Use a new device and pipette for each sample tested.
- Do not use past the expiration date of the serial.

### SAMPLES

- Acceptable sample types are feline serum, plasma, and anticoagulated (EDTA or Heparin) whole blood.
- Samples must be at room temperature and mixed gently before testing.
- Samples should be tested fresh or refrigerated at 2-8°C for up to one week.
- For longer storage, serum or plasma samples may be stored at -20°C. Frozen samples should be centrifuged before use.
- Avoid heavily hemolyzed or contaminated samples. Lipemic samples will not affect test results.

### INSTRUCTIONS FOR USE

- The kit and all samples to be tested should be brought to room temperature for approximately 30 minutes prior to testing.
- Remove the test device from the protective pouch and place horizontally on a flat surface.
- Using the provided pipette, vertically dispense **1 drop of sample** into each of the sample wells without touching the pipette to the test device. Allow sample to absorb.

 **X1 Each**



- Add **2 drops of chase buffer** vertically to the sample window.


 **X2 Each**




- Interpret the results at **10 minutes**. 

### RESULTS INTERPRETATION


- The test is positive if two lines are present (a control line and a test line) on the same test strip (FIV Antibodies or FeLV antigen positive).




**FIV Positive**



**FeLV Positive**




**Both Tests Positive**











**Both Tests Negative**

- The test is invalid if no line appears at the control line even if a line is present at the test line.



- All results should be interpreted within the context of all available clinical and diagnostic information.

### SYMBOL KEY

	Use by date		Temperature Limitation
	Batch code (lot)		Consult Instructions for use
	Serial number		Manufacturer
	Catalog number		Sufficient for <n> tests